

## Product Datasheet

### Dinitrophenyl hapten Antibody (Biotin) (orb1735923)

|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Description</b>         | Dinitrophenyl hapten Antibody (Biotin)                |
| <b>Species/Host</b>        | Rat                                                   |
| <b>Conjugation</b>         | Biotin                                                |
| <b>Tested Applications</b> | ELISA, WB                                             |
| <b>Target</b>              | Dinitrophenyl hapten                                  |
| <b>Preservatives</b>       | PBS/0,1%NaN3/50%glycérol.<br>PBS/0,1%NaN3/50%glycérol |
| <b>Form/Appearance</b>     | Biotin                                                |
| <b>Concentration</b>       | 1mg/mL                                                |
| <b>Storage</b>             | -20°C                                                 |
| <b>Note</b>                | For research use only                                 |
| <b>Isotype</b>             | IgM kappa                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                       |