

Product Datasheet

Dinitrophenyl hapten Antibody (orb1735935)

Description	Dinitrophenyl hapten Antibody
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Dinitrophenyl hapten
Preservatives	PBS/0,1%NaN ₃ . PBS/0,1%NaN ₃
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Expiration Date	12 months from date of receipt.