

Product Datasheet

Dinitrophenyl hapten Antibody (FITC) (orb1735981)

Description	Dinitrophenyl hapten Antibody (FITC)
Species/Host	Rat
Conjugation	FITC
Tested Applications	ELISA, WB
Target	Dinitrophenyl hapten
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	FITC
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Expiration Date	12 months from date of receipt.