

Product Datasheet

Rat IgG2c kappa (Isotype control) Antibody (orb1736031)

Description	Rat IgG2c kappa (Isotype control) Antibody
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IA
Preservatives	PBS/0,1%NaN ₃ . PBS/0,1%NaN ₃
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2c kappa
Expiration Date	12 months from date of receipt.