

Product Datasheet

Camelid conserved sequence Antibody (Biotin) (orb1736354)

Description	Rabbit Anti-Camelid VHH Antibody [Biotin], mAb is specific for Camelid VHH Antibody (Alpaca, Llama a...
Conjugation	Unconjugated
Tested Applications	FC, IF
Target	Camelid conserved sequence
Preservatives	33 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	- 20. Shelf life: 1 year
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Expiration Date	12 months from date of receipt.