

Product Datasheet

Camelid conserved sequence Antibody (orb1736355)

Description	Rabbit Anti-Camelid VHH Antibody, mAb is specific for Camelid sdAb (Alpaca, Llama and Camel). The mo...
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Camelid conserved sequence
Preservatives	32 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	- 20. Shelf life: 1 year
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Expiration Date	12 months from date of receipt.