

Product Datasheet

Camelid conserved sequence Antibody (orb1736357)

Description	Rabbit Anti-Camelid VHH Cocktail [iFluor 555] is specific for Camelid VHH antibody (Llama, Alpaca an...
Conjugation	Unconjugated
Tested Applications	FC, IF
Target	Camelid conserved sequence
Preservatives	30 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	- 20. Shelf life: 1 year
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Expiration Date	12 months from date of receipt.