

Product Datasheet

GFP and other derivatives of GFP Antibody (orb1736385)

Description	Primary vNAR single domain antibody
Species/Host	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GFP and other derivatives of GFP
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin
Form/Appearance	Liquid
Concentration	1mg/mL
Storage	Storage Condition: 4°C; Do not freeze. Shelf Life: 1 year
Note	For research use only
Isotype	Shark vNAR
Clonality	Monoclonal
Purity	> 95% as determined by SDS-PAGE
MW	24 kDa
Expiration Date	12 months from date of receipt.