

## Product Datasheet

### anti-PCSK9 VHH antibody (orb1736397)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Primary VHH single domain antibody                           |
| <b>Species/Host</b>        | Camelus                                                      |
| <b>Conjugation</b>         | Unconjugated                                                 |
| <b>Tested Applications</b> | IF, IHC                                                      |
| <b>Target</b>              | Human Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9)  |
| <b>Preservatives</b>       | 25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin    |
| <b>Form/Appearance</b>     | Liquid                                                       |
| <b>Concentration</b>       | 1mg/mL                                                       |
| <b>Storage</b>             | Storage Condition: 4°C; Do not freeze.<br>Shelf Life: 1 year |
| <b>Note</b>                | For research use only                                        |
| <b>Isotype</b>             | Camelid VHH                                                  |
| <b>Clonality</b>           | Monoclonal                                                   |
| <b>Purity</b>              | > 95% as determined by SDS-PAGE                              |
| <b>MW</b>                  | 15 kDa                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                              |