

Product Datasheet

Mouse IL1b(Interleukin 1 Beta) Microsample ELISA Kit (orb1736445)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P10749
Sensitivity	6.4 pg/mL
Expiration Date	Please enquire.

