

Product Datasheet

Rat PIK3C2a(Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit (orb1736455)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Tissue homogenates and other biological fluids.
Uniprot ID	A0A0G2K344
Sensitivity	0.071 ng/mL
Expiration Date	Please enquire.

