

## Product Datasheet

### Goat Cor(Cortisol) ELISA Kit (orb1736492)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| <b>Reactivity</b>      | Goat                                                                                                    |
| <b>Range</b>           | 3.13-200 ng/mL                                                                                          |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Sample Types</b>    | serum, plasma, urine and other biological fluids                                                        |
| <b>Sensitivity</b>     | 0.93 ng/mL                                                                                              |
| <b>Expiration Date</b> | Please enquire.                                                                                         |

