

Product Datasheet

Mouse SSC4D(Scavenger Receptor Cysteine Rich Family Member with 4 Domains) ELISA Kit (orb1736503)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	A1L0T3
Sensitivity	0.063 ng/mL
Expiration Date	Please enquire.

