

Product Datasheet

Human PHGDH(Phosphoglycerate Dehydrogenase) ELISA Kit (orb1736510)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.32-20 ng/mL
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	O43175
Sensitivity	0.14 ng/mL
Expiration Date	Please enquire.

