

Product Datasheet

Human 5-MTHF(5-Methyltetrahydrofolate) ELISA Kit (orb1736511)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P49914
Sensitivity	0.061 ng/mL
Expiration Date	Please enquire.

