

Product Datasheet

Human sIL2Ra(Soluble Interleukin-2 Receptor alpha) ELISA Kit (orb1736520)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	31.25-2000 pg/mL
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P01589
Sensitivity	17.1 pg/mL
Expiration Date	Please enquire.

