

Product Datasheet

Human CERS2(Ceramide synthase 2) ELISA Kit (orb1736521)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q96G23
Sensitivity	0.051 ng/mL
Expiration Date	Please enquire.

