

Product Datasheet

Human TNFaIP6(Tumor Necrosis Factor Alpha Induced Protein 6) ELISA Kit (orb1736522)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	P98066
Sensitivity	0.053 ng/mL
Expiration Date	Please enquire.

