

Product Datasheet

Rat ESM1(Endothelial Cell Specific Molecule 1) ELISA Kit (orb1736524)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	78.13-5000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P97682
Sensitivity	6.2 pg/mL
Expiration Date	Please enquire.

