

Product Datasheet

Rat TRPV4(Transient Receptor Potential Cation Channel Subfamily V, Member 4) ELISA Kit (orb1736525)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	46.88-3000 pg/mL
Note	For research use only
Sample Types	Tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9ERZ8
Sensitivity	18.87 pg/mL
Expiration Date	Please enquire.

