

Product Datasheet

Rat NBC4(Electrogenic Sodium Bicarbonate Cotransporter 4) ELISA Kit (orb1736542)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.79-50 ng/mL
Note	For research use only
Sample Types	Tissue homogenates and other biological fluids.
Uniprot ID	Q6RI88
Sensitivity	0.31 ng/mL
Expiration Date	Please enquire.

