

Product Datasheet

Rat ERFE(Erythroferrone) ELISA Kit (orb1736562)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	D4AB34
Sensitivity	0.071 ng/mL
Expiration Date	Please enquire.

