

Product Datasheet

Chicken GSTs(Glutathione S Transferases) ELISA Kit (orb1736614)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Gallus
Range	1.57-100 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O73888
Sensitivity	0.61 ng/mL
Expiration Date	Please enquire.

