

Product Datasheet

Chicken PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit (orb1736615)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Gallus
Range	3.13-200 ng/mL
Note	For research use only
Sample Types	tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9DEA3
Sensitivity	1.31 ng/mL
Expiration Date	Please enquire.

