

Product Datasheet

Human Rock-2(Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit (orb1736623)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.32-20 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	075116
Sensitivity	0.14 ng/mL
Expiration Date	Please enquire.

