

Product Datasheet

Human GCC2(Grip And Coiled-Coil Domain-Containing Protein 2) ELISA Kit (orb1736630)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	Q8IWJ2
Sensitivity	0.061 ng/mL
Expiration Date	Please enquire.

