

Product Datasheet

Human CART(Cocaine And Amphetamine Regulated Transcript) ELISA Kit [Out of stock] (orb1736641)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	31.25-2000 pg/mL
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	Q16568
Sensitivity	18.6 pg/mL
Expiration Date	Please enquire.

