

Product Datasheet

Mouse YM1(Chitinase 3-like 3) ELISA Kit (orb1736685)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	62.5-4000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O35744
Sensitivity	38.1 pg/mL
Expiration Date	Please enquire.

