

## Product Datasheet

### Chicken CCK(Cholecystokinin) ELISA Kit (orb1736686)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| <b>Reactivity</b>      | Gallus                                                                                                  |
| <b>Range</b>           | 15.63-1000 pg/mL                                                                                        |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Sample Types</b>    | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids    |
| <b>Uniprot ID</b>      | <b>Q9PU41</b>                                                                                           |
| <b>Sensitivity</b>     | 4.82 pg/mL                                                                                              |
| <b>Expiration Date</b> | Please enquire.                                                                                         |

