

Product Datasheet

Human NEU3(Sialidase-3) ELISA Kit (orb1736689)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.24-15 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9UQ49
Sensitivity	0.091 ng/mL
Expiration Date	Please enquire.

