

Product Datasheet

Human NEU4(Sialidase-4) ELISA Kit (orb1736690)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	93.75-6000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8WWR8
Sensitivity	39.1 pg/mL
Expiration Date	Please enquire.

