

## Product Datasheet

### Mouse GAD2(Glutamate Decarboxylase 2, Acid) ELISA Kit (orb1736721)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>      | Mouse                                                                                                   |
| <b>Range</b>           | 78.13-5000 pg/mL                                                                                        |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Sample Types</b>    | tissue homogenates, cell lysates, cell culture supernates and other biological fluids                   |
| <b>Sensitivity</b>     | 31.2 pg/mL                                                                                              |
| <b>Expiration Date</b> | Please enquire.                                                                                         |

