

Product Datasheet

Mouse SLC6A1(Sodium- and chloride-dependent GABA transporter 1) ELISA Kit (orb1736722)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P31648
Sensitivity	0.074 ng/mL
Expiration Date	Please enquire.

