

Product Datasheet

Human SARM1(Sterile Alpha And TIR Motif-Containing Protein 1) ELISA Kit (orb1736731)

Description The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

Reactivity Human

Range 0.16-10 ng/mL

Note For research use only

Sample Types serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Uniprot ID [Q6SZW1](#)

Sensitivity 0.071 ng/mL

Expiration Date Please enquire.

