

Product Datasheet

Mouse AHRR(Aryl Hydrocarbon Receptor Repressor) ELISA Kit (orb1736733)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Tissue homogenates, cell lysates and other biological fluids
Uniprot ID	Q3U1U7
Sensitivity	0.097 ng/mL
Expiration Date	Please enquire.

