

Product Datasheet

Human EIF2A(Eukaryotic Translation Initiation Factor 2 subunit alpha) ELISA Kit (orb1736765)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	78.13-5000 pg/mL
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids
Uniprot ID	Q9BY44
Sensitivity	31.5 pg/mL
Expiration Date	Please enquire.

