

Product Datasheet

SARS-CoV-2 Spike Protein S1 Recombinant Antigen (His tag) Ultra (orb1736821)

Description	Recombinant SARS-CoV-2 Spike protein (S1) fused to a 6X Histidine tag. Produced in Nicotiana bentham...
Species/Host	Plant
Tested Applications	ELISA, WB
Preservatives	0.1 M Phosphate Buffered Saline (PBS), pH = 7.4
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	-20°C Freezer
Note	For research use only
Formula	Monomer
Isotype	SARS-CoV-2
Clonality	Monoclonal
Purity	> 90 % as determined by SDS-PAGE.
Source	SARS-CoV-2
MW	± 123 kDa
Uniprot ID	P59594
Hazard Information	Non-Hazardous
Dilution Range	Optimal dilutions/concentrations should be determined by the user.
Expiration Date	6 months from date of receipt.