

Product Datasheet

PtX Mouse Anti-cytotoxic T-lymphocyte-associated protein 4 Recombinant Antibody (orb1736837)

Description

Recombinant mouse monoclonal

	lymphocyte-associated protein 4 (CTLA-4). ...
Species/Host	Plant
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Cytotoxic T-lymphocyte-associated protein 4 (CTLA-4)
Preservatives	0.1 M Phosphate Buffered Saline (PBS), pH = 7.4
Form/Appearance	Liquid
Concentration	1,0 mg/ml
Storage	-20°C Freezer
Note	For research use only
Formula	Full-length IgG
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 % as determined by SDS-PAGE
Source	Mouse
MW	± 150 kDa
Uniprot ID	P16410
Sensitivity	Refer to ELISA dose response graph in Datasheet for ELISA sensitivity.
Hazard Information	Non-Hazardous
Dilution Range	Suggested dilution is 1: 2 000-1: 10 000 for ELISA. Optimal dilutions/concentrations should be determined by the user.
Expiration Date	12 months from date of receipt.