

Product Datasheet

**SARS-CoV-2 Spike protein (RBD of S1 subunit)
Antibody (orb1736854)**

Description

Recombinant human monoclonal

Binding Domain of SARS-CoV-2 Spike protein...

Species/Host Plant

Conjugation Unconjugated

Tested Applications ELISA, WB

Target SARS-CoV-2 Spike protein (RBD of S1 subunit)

Preservatives 0.1 M Phosphate Buffered Saline (PBS), pH = 7.4

Form/Appearance Liquid

Concentration 1.0 mg/ml

Storage -20°C Freezer

Note For research use only

Formula Full-length IgG

Isotype IgG1

Clonality Monoclonal

Purity 99.30 % as determined by Mass Spectrometry.

Source Human

MW ± 150 kDa

Uniprot ID [P59594](#)

Hazard Information Non-Hazardous

Dilution Range Optimal dilutions/concentrations should be determined by the user.

Expiration Date 12 months from date of receipt.