

Product Datasheet

angiopoietin like 3 Antibody (orb1741746)

Description

angiopoietin like 3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human ANGPTL3. The exact sequence is proprietary.

Target

angiopoietin like 3

Preservatives

No Preservative. PBS

Form/Appearance

Liquid

Concentration

1.4 mg/ml (Please refer to the vial label for the specific concentration.)

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

54

Uniprot ID

[Q9Y5C1](#)

Dilution Range

WB: 1:500-1:3000. IHC-P: 1:100-1:1000.
*Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.

Expiration Date

12 months from date of receipt.