

Product Datasheet

neuropilin 1 Antibody (orb1742093)

Custom Reagents.
Description

neuropilin 1 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human Neuropilin 1. The exact sequence is proprietary.
Target	neuropilin 1
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol
Form/Appearance	Liquid
Concentration	0.79 mg/ml (Please refer to the vial label for the specific concentration.)
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: HeLa
Isotype	IgG
Clonality	Polyclonal
MW	103
Uniprot ID	O14786
Dilution Range	WB: 1:500-1:3000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.
Expiration Date	12 months from date of receipt.