

Product Datasheet

ribosomal protein lateral stalk subunit P2

Antibody (orb1742278)

Description

ribosomal protein lateral stalk subunit

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human RPLP2. The exact sequence is proprietary.
Target	ribosomal protein lateral stalk subunit P2
Preservatives	0.01% Thimerosal. 1XPBS (pH7), 1% rAlbumin, 20% Glycerol
Form/Appearance	Liquid
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: 293T
Isotype	IgG
Clonality	Polyclonal
MW	12
Uniprot ID	P05387
Dilution Range	WB: 1:5000-1:10000. ICC/IF: 1:100-1:1000. IHC-P: 1:100-1:1000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.
Expiration Date	12 months from date of receipt.