

Product Datasheet

ATPase H⁺ transporting V1 subunit E1 Antibody (orb1742290)

Description

ATPase H⁺ transporting V1 subunit E1

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human ATP6V1E1. The exact sequence is proprietary.
Target	ATPase H ⁺ transporting V1 subunit E1
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol
Form/Appearance	Liquid
Concentration	0.17 mg/ml (Please refer to the vial label for the specific concentration.)
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: Mouse brain , HepG2 , Molt-4
Isotype	IgG
Clonality	Polyclonal
MW	26
Uniprot ID	P36543
Dilution Range	WB: 1:500-1:3000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.
Expiration Date	12 months from date of receipt.