

Product Datasheet

poly(A) binding protein nuclear 1 Antibody (orb1742380)

Description

poly(A) binding protein nuclear 1

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human PABPN1. The exact sequence is proprietary.
Target	poly(A) binding protein nuclear 1
Preservatives	0.01% Thimerosal. 1XPBS (pH7), 1% rAlbumin, 20% Glycerol
Form/Appearance	Liquid
Concentration	0.49 mg/ml (Please refer to the vial label for the specific concentration.)
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: mouse brain , Jurkat
Isotype	IgG
Clonality	Polyclonal
MW	33
Uniprot ID	Q86U42
Dilution Range	WB: 1:500-1:3000. ICC/IF: 1:100-1:1000. IHC-P: 1:100-1:1000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.
Expiration Date	12 months from date of receipt.