

Product Datasheet

Rho GDP dissociation inhibitor alpha Antibody (orb1742404)

Description

Rho GDP dissociation inhibitor alpha

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human Rho GDI alpha. The exact sequence is proprietary.
Target	Rho GDP dissociation inhibitor alpha
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol
Form/Appearance	Liquid
Concentration	0.76 mg/ml (Please refer to the vial label for the specific concentration.)
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: Mouse brain
Isotype	IgG
Clonality	Polyclonal
MW	23
Uniprot ID	P52565
Dilution Range	WB: 1:500-1:10000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.
Expiration Date	12 months from date of receipt.