

Product Datasheet

forkhead box M1 Antibody (orb1742698)

Description

forkhead box M1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ChIP, WB

Applications
Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human FOXM1. The exact sequence is proprietary.

Target

forkhead box M1

Preservatives

0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid

Concentration

1.52 mg/ml (Please refer to the vial label for the specific concentration.)

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: FOXM1-transfected 293T

Isotype

IgG

Clonality

Polyclonal

MW

84

Uniprot ID
[Q08050](https://www.uniprot.org/entry/Q08050)
Dilution Range

WB: 1:1000-1:10000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.

Expiration Date

12 months from date of receipt.