

Product Datasheet

forkhead box D1 Antibody (orb1742700)

Description

forkhead box D1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human FOXD1. The exact sequence is proprietary.

Target

forkhead box D1

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid

Concentration

1 mg/ml (Please refer to the vial label for the specific concentration.)

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: 293T nuclear extract

Isotype

IgG

Clonality

Polyclonal

MW

46

Uniprot ID
[Q16676](#)
Dilution Range

WB: 1:500-1:3000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.

Expiration Date

12 months from date of receipt.