

Product Datasheet

**FGR proto-oncogene, Src family tyrosine kinase
Antibody (orb1742702)**

Description

FGR proto-oncogene, Src family

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human FGR. The exact sequence is proprietary.
Target	FGR proto-oncogene, Src family tyrosine kinase
Preservatives	0.01% Thimerosal. 1XPBS (pH7), 20% Glycerol
Form/Appearance	Liquid
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: HeLa , HepG2 , Molt-4
Isotype	IgG
Clonality	Polyclonal
MW	59
Uniprot ID	P09769
Dilution Range	WB: 1:500-1:3000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.
Expiration Date	12 months from date of receipt.