

Product Datasheet

Fc fragment of IgG receptor IIb Antibody (orb1742703)

Description

Fc fragment of IgG receptor IIb Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, ICC, WB

Applications
Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human CD32B. The exact sequence is proprietary.

Target

Fc fragment of IgG receptor IIb

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid

Concentration

1 mg/ml (Please refer to the vial label for the specific concentration.)

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: THP-1

Isotype

IgG

Clonality

Polyclonal

MW

34

Uniprot ID

P31994

Dilution Range

WB: 1:500-1:3000. ICC/IF: 1:100-1:1000.
*Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.

Expiration Date

12 months from date of receipt.