

Product Datasheet

**nuclear receptor subfamily 2 group F member 6
Antibody (orb1742705)**

Description

nuclear receptor subfamily 2 group F

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ChIP, IHC-P, WB

Applications
Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human NR2F6. The exact sequence is proprietary.

Target

nuclear receptor subfamily 2 group F member 6

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid

Concentration

1 mg/ml (Please refer to the vial label for the specific concentration.)

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

43

Uniprot ID
P10588
Dilution Range

WB: 1:500-1:3000. IHC-P: 1:100-1:1000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.

Expiration Date

12 months from date of receipt.